

Display Report

Analysis Info

Analysis Name D:\Data\test\0115\qiangjilaochun_56_01_2752.d
Method 2752.m
Sample Name qiangjilaochun
Comment

Acquisition Date 1/15/2016 11:58:50 AM

Operator brucker
Instrument amaZon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Positive	Alternating Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	1000 m/z
Capillary Exit	140.0 V	n/a	n/a	Trap Drive	50.6
Accumulation Time	2257 μ s	Averages	5 Spectra	Auto MS/MS	off

